

GENERAL CHEMISTRY LAB COMPONENT CHE101L CONTENT: LAB 3

Dissolution Reactions: Heats of Dissociation

NAME

Joy kuman Ghosh

SECTION 7

STUDENT ID

DATE & TIME .

22 11 424 642

21.12.2023 11:15 am

NAME OF THE INSTRUCTOR

Justo)

SIGNATURE & DATE

07.01.2024

REPORT SUBMISSION DATE

EXPERIMENT 3

SESSION 1 (QUALITATIVE)

DISSOLUTION REACTIONS: HEATS OF DISSOCIATION

Heats (exothermic or endothermic) are associated with chemical reactions. Quantity of heat evolved or absorbed is directly proportional to the amount reacted. Consider the reactions below:

MgSO4(s) + H₂O (I)
Mg
$$^{2+}(aq.) + SO_4 ^{2-}(aq.) \pm heat$$

NaNO₃(s) + H₂O (I)
Na⁺(aq.) + NO₃⁻(aq.) $\pm heat$

Heat could be generated or absorbed in this reaction. When heat is generated/released from a chemical reaction it is called exothermic reaction (you can feel it by touching the reaction container (warmer) and when heat is absorbed the reaction is called endothermic (colder). When reactions occur in a reaction vessel (e.g., Beaker) in aqueous condition, formation and dissociation of chemical bonds occur simultaneously. Bond formation and dissociation involves heat energy of the system which is expressed by the term Q which is called enthalpy.

PROBLEM STATEMENT: Is heat energy related to chemical reactions, how?

This experiment is subdivided into two parts:

- I. QUALITATIVE
- II. QUANTITATIVE

PART I. QUALITATIVE

DATA COLLECTION:

Place about 30 mL of distilled water into a 50 mL beaker. Suspend a thermometer (having 0.1°C division mark) into the beaker using thermometer clamp and ring stand. Please make sure that the thermometer is not touching the bottom of the beaker, as any movement of the beaker could break the thermometer. Record the temperature of water in the beaker in every 30 seconds for 180 seconds.




FIGURE 1: Experimental setup for dissolution reaction

Place a moderate amount (which would be 1 to 3 cm³) of supplied anhydrous magnesium sulfate (MgSO₄) to the beaker. Mix vigorously with the glass rod for 2 minutes. Record your observations. (2 points)

We observed distilled water stability for 60 seconds, and the temperature of the water was 24'10°C, then we added MgSOy. After 120 seconds, we noted the temperature increase to 34°C, indicating an exothermic neaction.

Repeat this procedure with each of the following compounds with two different amounts (roughly 1:2): (2 points)

- a. Sodium Nitrate, NaNO3 24.25°C 23°C Endothermic
- b. Sodium Chloride, NaCl 24°C 23°C Endo thermic
- c. Hydrated Calcium Chloride, CaCl2.2H2O -24°C 26°C Enothermic
- d. Ammonium Nitrate, NH4NO3 24°C 22°C Endo thermic
- e. Magnesium sulfate, MgSO4 24'10°C 34°C Enothermic

DATA ANALYSIS:

Compare and contrast the behaviors of these compounds. Identify any generalizations that can be made about the chemical reactions observed. What conclusions can be drawn from the data? (4 points)

The compounds that we've used in the experiment are NaNO3, NaCl, CaCl2. 2H2O, NH4NO3, MgSO4. Among these compounds, all of them are soluble in water.

Herre, NaNO3, Nacl, NH4NO3 absorb energy and decrease the temperature when dissolved in water. So, NaNO3, Nacl and NH4NO3 are Endothermic.

On the other hand, Mysoy and Cacl2. 2 Hzo Hearelease heat and rise the temperature when dissolved in water. Thus Mysoy and Cacl2. 2 HzO are Exothermic.

SESSION 2(QUANTITATIVE)

PART II. QUANTITATIVE

DATA COLLECTION:

a. Accurately weigh a 3 to 5 gm sample of MgSO₄ on the analytical balance. Record the exact mass here. For 4 different trials below measure four different weight samples (e.g., 1,2, 4 & 5 grams respectively).

b. Suspend the thermometer into a polystyrene cup/coffee cup. Make sure of the thermometer is not touching the bottom of the cup. Measure 20 mL of distilled water by a volumetric cylinder into the cup and stir for 240 second. Record the temperature in every 20 seconds. After 240 seconds add MgSO₄ with vigorous mixing while continuing to record data for 5 minutes.

c. Determine the temperature change, ΔT , for the reaction. This can be done from the difference of the highest temperature minus the slope of the line go through the points from first 240 seconds of data.

d. Draw a temperature vs. time graph. Draw the best curve through the points and point out what is happing in each part of the curve.

e. Repeat the whole procedure with NaNO_{3.}

DATA TABLE 1:

TRAILS

(I) Mass of MgSO ₄ _0.5 gm		(II) Mass of N	(II) Mass of MgSO ₄ _1.0 gm	
Time (s)	Temp(°C)	Time(s)	Temp(°C)	
20	25	20	25	
40	25	40	25.10	
60	25	60	25.5	
80	27	80	29	
100	27	100	31	
120	26.35	120	31.10	
140	26.5	140	31.10	
160	26.25	160	31	
180	26	180	34	
200	26	200	31	
220	25.9	220	31	
240	25.75	240	30.9	
260	25.5	260	30.5	
280	25.5	280	30.5	
300	25.5	300	30.25	

(III) Mass of MgSO ₄ _1 <u>.5 gm</u>		(IV) Mass of I	(IV) Mass of MgSO _{4_2.0 gm}	
Time (s)	Temp(°C)	Time(s)	Temp(°C)	
20	23	20	24.5	
40	23	40	24.5	
60	23	60	24.5	
80	29.5	80	30	
100	32	100	33	
120	32	120	33.5	
140	32	140	34	
160	32	160	34	
180	31:5	180	34	
200	31:5	200	33.8	
220	31.5	220	33.8	
240	31	240	33.8	
260	31	260	33.5	
280	31	280	33.5	
300	30.75	300	33.25	

DATA TABLE 2:

TRAILS

(I) Mass of NaNO ₃ _0.5 gm		(II) Mass of N	(II) Mass of NaNO _{3_} 1.00 gm	
Time (s)	Temp(°C)	Time(s)	Temp(°C)	
20	24.75	20	23	
40	24.75	40	23	
60	24.75	60	23	
80	24.5	80	22:5	
100	24.25	100	22	
120	24.25	120	22	
140	24.1	140	22	
160	24.1	160	22	
180	24.1	180	22	
200	24	200	21.75	
220	24	220	21.35	
240	24	240	21.75	
260	24	260	21.5	
280	24	280	21:5	
300	24.1	300	21.5	

(III) Mass of NaNO ₃ _1.5 gm		(IV) Mass of I	(IV) Mass of NaNO ₃ _2.0 gm	
Time (s)	Temp(°C)	Time(s)	Temp(°C)	
20	25.25	20	24.25	
40	25.10	40	24.1	
60	25.10	60	24	
80	25	80	23.75	
100	24.5	100	23	
120	24'1	120	23	
140	241	140	23	
160	24	160	23	
180	24	180	23	
200	23'9	200	23.25	
220	23.9	220	23.5	
240	23.9	240	22	
260	23.75	260	22:5	
280	23.75	280	22	
300	23.75	300	21.2	

DATA ANALYSIS

1. What do you understand form both data sets you recorded and from the other trials? (4 points)

For MgSOy: From the data tables of MgSOy, we can say that MgSoy dissolves in water & increase temperature for certain time period and then it stanted to decnease. Thus it's an Exothermic reaction.

 $MgSO_4(s) + H_2O(L) \rightarrow Mg^{2+}(a_2) + SO_4^{2-}(a_2)$, $\Delta H = -Ve$ (Fnothermic

If we add mone substance then we can see the mone temperature to TUSE.

For NaNOz: From the data tables of NaNOz, we can say that NaNOz dissolves in water and decrease temperature for a centain time period and then it stanted increasing again. This it's an Endothermic.

 $NaN0_3(s) + H_20(\lambda) \longrightarrow Na^{\dagger}(a_2) + No_3^{\dagger}(a_2)$, $\Delta = +ve$ (Endo-thenmic Reaction)

If we add mone substance, then we can see the mone tempenature to fall.

2. Calculate the heat, Q & moles, n, of the reaction both data sets. Take help from the equitation Q = C x M x Δ T. Assume C = 4.18 Joules/gram 0 C and M is the mass of water (take the water density as 1.00 grams/cm³). (4 points)

For MgSO4!

0'5 gm: Heat,
$$Q = C \times M \times \Delta T = 4'18 \times 20 \times (27-25) = 167.2 \text{ j}$$

Mole, $n = \frac{0.5}{120} = 0.0041 \text{ mole}$

1 gm: Heat, $Q = C \times M \times \Delta T = 4'18 \times 20 \times (31'10-25.5) = 468'16 \text{ j}$

Mole, $n = \frac{1}{120} = 0.0093 \text{ mole}$

1'5 gm: Heat, $Q = C \times M \times \Delta T = 4'18 \times 20 \times (31'10-25.5) = 468'16 \text{ j}$

Mole, $n = \frac{1.5}{120} = 0.0093 \text{ mole}$

1'5 gm: Heat, $Q = C \times M \times \Delta T = 4'18 \times 20 \times (32-23) = 752'4 \text{ j}$

Mole, $n = \frac{1.5}{120} = 0.0125 \text{ mole}$

2gm: Heat, $Q = C \times M \times \Delta T = 4'18 \times 20 \times (34-24'5) = 794'2 \text{ j}$

Mole, $n = \frac{2}{120} = 0.0167 \text{ mole}$

1 gm: Heat, $Q = C \times M \times \Delta T = 4'18 \times 20 \times (24'75-24) = 62'7 \text{ j}$

Mole, $n = \frac{0.5}{35} = 0.0059 \text{ mole}$

1 gm: Heat, $Q = C \times M \times \Delta T = 4'18 \times 20 \times (23'10-23'75) = 112'8' \text{ j}$

Mole, $n = \frac{1}{35} = 0.0118 \text{ mole}$

1 5gm: Heat, $Q = C \times M \times \Delta T = 4'18 \times 20 \times (25'10-23'75) = 112'8' \text{ j}$

Mole, $n = \frac{1.5}{35} = 0.0176 \text{ mole}$


2 gm: Heat, $Q = C \times M \times \Delta T = 4'18 \times 20 \times (24'25) = 209 \text{ j}$

Mole, $n = \frac{1.5}{35} = 0.0176 \text{ mole}$

2 gm: Heat, $Q = C \times M \times \Delta T = 4'18 \times 20 \times (24'25) = 209 \text{ j}$

Mole, $n = \frac{1.5}{35} = 0.0176 \text{ mole}$

3. Plot the collected data as moles, n vs. Q both sets of data. Number of moles can be calculated as n = (mass of sample in gram) / (molecular weight in grams/mole). Try to



MENTAL MODEL: Use the chemical equation given above to represent the dissolution reaction in this experiment. Draw a picture(s) which describes what is happening in atomic or in molecular level for either MgSO₄ or NaNO₃ system. How heat release or absorbed can be described from these pictures? (5 points)

The formula for MgSO4 dissolution: $MgSO4(s) + H_{2}O(N) \longrightarrow Mg^{2+}(aq) + SO_{4}^{2}(aq) + H_{2}O(A)$

The dissolution of MgSo4 in water is an exothermic process, meaning that it releases heat. This heat release is due to the strong affiniactions between the ions in MgSo4 and the water molecules. When MgSo4 is dissolved in water, the ions dissociate and break apart, and water molecules surround and stabilize them. This process of hydration releases a significant amount of energy in form of heat.